



Sequencing by binding (SBB[®]) delivers unprecedented NGS accuracy

June 8, 2022 | AGBT
Jonas Korlach | CSO

The last 18 months have been transformational for PacBio

Genomes

Epigenomes

Transcriptomes

Chromatin architecture & dynamics

Metagenomics

Gene therapy

SARS-CoV-2



Foundation for T2T, calling all variants

Simultaneous 5mC calling

MAS-Iso-Seq

Fiber-Seq/SAMOSA

Full-length 16S, complete MAGs

Complete AAV sequencing solution

HiFiViral

The last 18 months have been transformational for PacBio

Released new binding kits

Drove DNA input down >5×

Increased average yield >30%

Released new prep kits

Released automated protocols

Released HiFiViral kit

Enabled 5mC calling on instrument


Released AAV workflow

Enabled demultiplex on instrument

Empowered high-throughput processing

Consolidated workflows and protocols

And much more...



Accuracy matters — it's the hallmark of who we are

Human genetics — Neuroscience

Human genetics — Immunology

Rare + inherited disease research

Plant + animal sciences

Infectious disease / microbiology

Potential for early-stage cancer screening

Potential for cancer recurrence monitoring

Enabling therapy selection

Targeted clinical panels

Potential for noninvasive prenatal screening



HiFi sequencing


Delivers long reads with the highest accuracy — even in hard-to-sequence regions




SBB sequencing

Promises significant accuracy improvements over conventional NGS approaches

Going beyond what we shared at JP Morgan




Even lower error rate than original emulsion PCR method



Even lower error rate than original emulsion PCR method

Reaches 0.1% "true variation level" between the sample and reference



Sequencing by binding (SBB)

1

Technology

2

Benchmark

3

Application


Developing an innovative platform to house SBB

Novel and state-of-the-art inventions

4 core focus areas

140+ patents pending

~50 patents allowed/issued



Proprietary
chemistry



Consumables
design




Instrument
engineering



Software
algorithms


Sequencing around the clock

>1500 runs completed in 2022 alone



Breakthrough short-read sequencing

Key design principles and goals



Mid- to high-throughput NGS platform



Optical and mechanical innovations



Scalable, flexible, and cost-optimized



Unparalleled accuracy from SBB



SBB is fundamentally designed to maximize accuracy

SBB chemistry separates
interrogation and *incorporation* steps



Multiple optimization points
increase accuracy and flexibility

Blocked 3' end



Interrogate

Flow nucleotides, image, wash



Activate

Remove 3' RT



Incorporate


Flow blocked nucleotides




SBB advantages: Incorporates native nucleotides, produces unmodified DNA

No base modifications, no molecular scarring

Sequencing by synthesis (SBS)



Sequencing by binding (SBB)



Benefits of SBB over traditional short read NGS


>90% bases at Q40+

Low duplications rate

No index hopping

Sequence through difficult / repetitive regions

SBB offers near “perfect” reads




SBB offers best-in-class accuracy

SBB error rates ~15× lower at any given cycle, between 1:10,000 to 1:100,000

The complete sequence of a human genome

Sergey Nurk^{1,†}, Sergey Koren^{1,†}, Arang Rhie^{1,†}, Mikko Rautiainen^{1,†}, Andrey V. Bzikadze², Alla Mikheenko³, Mitchell R. Vollger⁴, Nicolas Altemose⁵, Lev Uralsky^{6,7}, Ariel Gershman⁸, Sergey Aganezov⁹, Savannah J. Hoyt¹⁰, Mark Diekhans¹¹, Glennis A. Logsdon⁴, Michael Alonge⁹, Stylianos E. Antonarakis¹², Matthew Borchers¹³, Gerard G. Bouffard¹⁴, Shellise Y. Brooks¹⁴, Gina V. Caldas¹⁵, Haoyu Cheng^{16,17}, Chen-Shan Chin¹⁸, William Chow¹⁹, Leonardo G. de Lima¹³, Philip C. Dishuck⁴, Richard Durbin²¹, Tatiana Dvorkina³, Ian T. Fiddes²², Giulio Formenti^{23,24}, Robert S. Fulton²⁵, Arkarachai Fungtammasan¹⁸, Erik Garrison^{11,26}, Patrick G.S. Grady¹⁰, Tina A. Graves-Lindsay²⁷, Ira M. Hall²⁸, Nancy F. Hansen²⁹, Gabrielle A. Hartley¹⁰, Marina Haukness¹¹, Kerstin Howe¹⁹, Michael W. Hunkapiller³⁰, Chirag Jain^{1,31}, Miten Jain¹¹, Erich D. Jarvis^{23,24}, Peter Kerpeldjiev³², Melanie Kirsche⁹, Mikhail Kolmogorov³³, Jonas Korlach³⁰, Milinn Kremitzki²⁷, Heng Li^{16,17}, Valerie V. Maduro³⁴, Tobias Marschall³⁵, Ann M. McCartney¹, Jennifer McDaniel³⁶, Danny E. Miller^{4,37}, James C. Mullikin^{14,29}, Eugene W. Myers³⁸, Nathan D. Olson³⁶, Benedict Paten¹¹, Paul Peluso³⁰, Pavel A. Pevzner³³, David Porubsky⁴, Tamara Potapova¹³, Evgeny I. Rogaev^{6,7,39,40}, Jeffrey A. Rosenfeld⁴¹, Steven L. Salzberg^{9,42}, Valerie A. Schneider⁴³, Fritz J. Sedlazeck⁴⁴, Kishwar Shafin¹¹, Colin J. Shew²⁰, Alaina Shumate⁴², Yumi Sims¹⁹, Arian F. A. Smit⁴⁵, Daniela C. Soto²⁰, Ivan Sović^{30,46}, Jessica M. Storer⁴⁵, Aaron Streets^{5,47}, Beth A. Sullivan⁴⁸, Françoise Thibaud-Nissen⁴³, James Torrance¹⁹, Justin Wagner³⁶, Brian P. Walenz¹, Aaron Wenger³⁰, Jonathan M. D. Wood¹⁹, Chunlin Xiao⁴³, Stephanie M. Yan⁴⁹, Alice C. Young¹⁴, Samantha Zarate⁹, Urvashi Surti⁵⁰, Rajiv C. McCoy¹⁹, Megan Y. Dennis²⁰, Ivan A. Alexandrov^{3,7,51}, Jennifer L. Gerton¹³, Rachel J. O'Neill¹⁰, Winston Timp^{8,42}, Justin M. Zook³⁶, Michael C. Schatz^{2,49}, Evan E. Eichler^{4,24,†}, Karen H. Miga^{11,†}, Adam M. Phillippy^{1,†}




SBB vs SBS empirical vs reported per-base Q score for CHM13

Uncalibrated SBB Q score correlates well with observed errors




Excellent variant calling performance for SBB



- Expecting further improvements through variant caller training

What does unprecedented accuracy look like?


Low-complexity region (28 As in the reference)



No high-quality reads
No variants called

What does unprecedented accuracy look like?

SBB cleanly sequences through 27 bp poly A (28 bp in the reference)



Example of a ‘difficult’ region – TOMM40



A TOMM40 variable-length polymorphism predicts the age of late-onset Alzheimer’s disease

AD Roses^{1,2}, MW Lutz^{1,2},
H Amrine-Madsen³,
AM Saunders^{1,2}, DG Crenshaw^{1,2},
SS Sundseth^{1,2}, MJ Huentelman⁴,
KA Welsh-Bohmer^{1,5} and
EM Reiman^{4,6,7}

¹Department of Medicine, Duke University,

The ε4 allele of the apolipoprotein E (APOE) gene is currently the strongest and most highly replicated genetic factor for risk and age of onset of late-onset Alzheimer’s disease (LOAD). Using phylogenetic analysis, we have identified a polymorphic poly-T variant, rs10524523, in the translocase of outer mitochondrial membrane 40 homolog (TOMM40) gene that provides greatly increased precision in the estimation of age of LOAD onset for APOE ε3 carriers. In two independent clinical cohorts, longer lengths of rs10524523 are associated with lower risk for LOAD in APOE ε3 carriers.

ARTICLE


TOMM40 ‘523’ poly-T repeat length is a determinant of longitudinal cognitive decline in Parkinson’s disease

Megan C. Bakeberg^{1,2}, Anastazja M. Gorecki^{1,3}, Abigail L. Pfaff^{1,4}, Madison E. Hoes¹, Sulev Köks^{1,4}, P. Anthony Akkari^{1,2,4}, Frank L. Mastaglia^{1,2,4} and Ryan S. Anderton^{1,2,5}

The translocase of outer mitochondrial membrane 40 (TOMM40 ‘523’ polymorphism has previously been associated with age of Alzheimer’s disease onset and cognitive functioning in non-pathological ageing, but has not been explored as a candidate risk marker for cognitive decline in Parkinson’s disease (PD). Therefore, this longitudinal study investigated the role of the ‘523’ variant in cognitive decline in a patient cohort from the Parkinson’s Progression Markers Initiative (PPMI). As such, a group of 360 people with PD were assessed annually for cognitive performance using multiple measures of cognition and neuroimaging. The results of the TOMM40 ‘523’ variant on cognitive performance were compared to the whole-genome data.


“Assay development for the TOMM40 ‘523’ variant is generally considered to be difficult, as poly-T variants are challenging to sequence.”¹

Example of a ‘difficult’ region – TOMM40




Example of a ‘difficult’ region – TOMM40

SBB cleanly sequences through poly-T repeat locus




Another similar region nearby





No high-quality reads
No variants called

Another similar region nearby




SBB correctly sequences through lengthy C/G repeats

Example Mucin 1 (*MUC1*) kidney disease



Sequencing performance on a healthy control sample


SBB did not observe spurious mutant variants in CHM13



Observed vs expected ctDNA at low variant allele frequency

Variant allele percent shows good linearity, even without use of UMIs

High sensitivity and specificity down to 0.125%, even with modest (<12,000×) coverage




Controls and reference materials
SeraSeq® ctDNA Mutation Mix v2

How low can SBB go without UMIs?

Tuberculosis amplicon shows near perfect linearity from 10% to 0.001%

	Res (C) counts	WT (G) counts	Total # counts	Observed % res
10%_rep1	2580549	18479388	21,059,937	12.2534%
10%_rep2	2998687	23062382	26,061,069	11.5064%
10%_rep3	743168	6508219	7,251,387	10.2486%
1%_rep1	76912	6388097	6,465,009	1.1897%
1%_rep2	46152	4025690	4,071,842	1.1334%
1%_rep3	65964	5334716	5,400,680	1.2214%
0.1%_rep1	22836	19283363	19,306,199	0.1183%
0.1%_rep2	8121	7646655	7,654,776	0.1061%
0.1%_rep3	6505	6587930	6,594,435	0.0986%
0.01%_rep1	1164	11078726	11,079,890	0.0105%
0.01%_rep2	197	2579126	2,579,323	0.0076%
0.01%_rep3	1030	10684737	10,685,767	0.0096%
0.001%_rep1	171	18047826	18,047,997	0.0009%
0.001%_rep2	117	8328474	8,328,591	0.0014%
0.001%_rep3	786	59804418	59,805,204	0.0013%



Where to from here for SBB?

Today

- Taking applications for collaboration; run your samples in our lab on SBB
- Visit our suite or www.pacb.com/sbb

Tomorrow

- Thurs, June 9, 8:00–8:30 am: *Advancing NGS accuracy by an order of magnitude*
Jennifer Stone, PhD, Vice President, Segment Marketing

Late Q3

- Formal external beta commences; more information to be shared at ASHG

1H 23

- On track for platform commercial availability

